

Thank you for helping the City of Los Angeles address supply needs for COVID-19.

The following packet provides the following information:

1. Manufacturing Specs
2. Patterns and Photos of Final Product
3. Box Labeling Requirements

Fabric Face Mask Material Specification

Fabric, Thread

Kaiser Permanente recognizes that fabric supplies are becoming increasingly tight, and we may need to accept alternatives. **However, all materials must be latex free.** Manufacturers should identify what they can source and communicate their options to Kaiser Permanente for approval. The general specification for mask materials includes:

FABRIC TYPE

- + Woven 6-8 oz (in densities similar to dress shirt or bed sheet fabrics)
- + Jersey Knit 6-8 oz
- + All masks will use two layers of fabric per the pattern and instructions

FABRIC COLOR

- + Lighter colors preferred

CONTENTS

- + 100% cotton or cotton blend

THREAD


- + Polyester per factory standards

COATINGS/TREATMENTS

- + No coatings or treatments required at this time including anti-microbial coatings

Large Mask

one piece of fabric 16 1/4 length and 8.5 width



Finished Masks



Packaging / Packing / Labeling of Bags

- + 504 masks per carton,
36 or 72 Polybags per carton
- + Package in standard approximate quart size
Polybag of 7 or 14 masks per bag with masks
stacked and folded
- + Sticker (ASE255) to be applied to outside
of each Polybag and to the carton with the
following text:

For Reusable Masks

FACE MASK
NOT INTENDED FOR CLINICAL USE
TO BE LAUNDERED BEFORE USE

For One-Time Use Masks

FACE MASK
NOT INTENDED FOR CLINICAL USE
NOT REUSABLE, DISCARD AFTER USE

Box Labeling Requirements

The following text must be visible on each box exterior. Select label option based on mask type.

For Reusable Masks

FACE MASK
NOT INTENDED FOR CLINICAL USE
TO BE LAUNDERED BEFORE USE

For One-Time Use Masks

FACE MASK
NOT INTENDED FOR CLINICAL USE
NOT REUSABLE, DISCARD AFTER USE